



MASKET CAPSULORHEXIS FORCEPS: OVERALL LENGTH 99 mm. CURVED SHAPED SHAFT LENGTH 13 mm. DELICATE CYSTOTOME SHAPED TIP 0.2 mm. (CL-FR-22)




UTRATA CAPSULORHEXIS FORCEPS: OVERALL LENGTH 99 mm. ANGLED STRAIGHT SHAFT LENGTH 13 mm. DELICATE CYSTOTOME SHAPED TIP 0.2 mm. (CL-FR-21)



JAFF TYING FORCEPS, STRAIGHT: OVERALL LENGTH 99 mm. 6 mm TYING PLATFORM. THREE HOLE BODY. DULL FINISH. (CL-FR-02)




JAFF TYING FORCEPS, CURVED: OVERALL LENGTH 99 mm. 6 mm TYING PLATFORM. THREE HOLE BODY. DULL FINISH. (CL-FR-03)



CASTROVIEJO TYING FORCEPS. STRAIGHT: OVERALL LENGTH 115 mm. THREE HOLE BODY. TYING PLATFORM 6 mm / 10 mm. DULL FINISH. (CL-FR-04)



CASTROVIEJO TYING FORCEPS. CURVED: OVERALL LENGTH 115 mm. THREE HOLE BODY. TYING PLATFORM 6 mm / 10 mm. DULL FINISH. SMOOTH CURVED. (CL-FR-05)




TYING FORCEPS CURVED. Sleek model: PERFORATED LIGHT WEIGHT MODEL. 8 mm TYING PLATFORM. ELEGANT LOOK. SMOOTH POLISHED 94 mm OVERALL LENGTH. (CL-FR-09)



TYING FORCEPS STRAIGHT. Sleek model: PERFORATED LIGHT WEIGHT MODEL ELEGANT LOOK. SMOOTH POLISHED 94 mm OVERALL LENGTH. 8 mm TYING PLATFORM. (CL-FR-09/S)




TYING FORCEPS, STRAIGHT : PERFORATED LIGHT WEIGHT MODEL. 87 mm OVERALL LENGTH.
6 mm TYING PLATFORM. (CL-FR-09/P)




BONN IRIS FORCEPS : OVERALL LENGTH 100 mm. 1X 2, 0.12 mm TEETH. THREE HOLE BODY.
(CL-FR-17)



KELMAN Mc PHERSON FORCEPS : OVERALL LENGTH 84 mm ANGLED SHAFT.
LONG 11 mm TYING PLATFORM. (CL-FR-24/H)



St.MARTIN FORCEPS: OVEALL LENGTH 85 mm. 1 X 2 TEETH. 6 mm TYING PLATFORM.
FLAT SERRATED HANDLE. STRAIGHT SHAFT. (CL-FR-12)



TISSUE FORCEPS, 1 X 2 TEETH: OVERALL LENGTH 100 mm. STRAIGHT SHAFT 1 X 2 TEETH. SERRATED HANDLE. (CL-FR-14)




BONN IRIS FORCEPS, Long handle: EXCLUSIVELY DESIGNED. OVERALL LENGTH 97 mm. 1 X 2, 0.12 mm TOOTH. BROAD 11 mm WIDE FLAT PERFORATED HANDLE FOR FIRM GRIP. (CL-FR-11/SP)




BONN IRIS FORCEPS EXCLUSIVELY DESIGNED: OVERALL LENGTH 97 mm. 1 X 2, 0.12 mm TOOTH. BROAD 11 mm WIDE FLAT PERFORATED HANDLE FOR FIRM GRIP. (CL-FR-11/ED)



COLIBRI FORCEPS, Perforated body: OVERALL LENGTH 78 mm. 1 X 2, 0.12 mm TEETH SET AT 45°. FLAT PERFORATED HANDLE. (CL-FR-19)




COLIBRI FORCEPS. *Serrated body:* OVERALL LENGTH 78 mm. 1X 2, 0.12 mm TEETH SET AT 45°. FLAT SERRATED HANDLE. (CL-FR-19/SR)



LIMS FORCEPS. OVERALL LENGTH 85 mm. 1X 2, 0.12 mm TEETH. FLAT SERRATED HANDLE. (CL-FR-18)




LENS HOLDING FORCEPS: FOR HOLDING & INSERTING FOLDABLE IOL. OVERALL LENGTH 118 mm. SERRATED, PERFORATED POLISHED BODY. (CL-FR-25)



ARRUGA CAPSULE FORCEPS : OVERALL LENGTH 100 mm. 2 mm CUP. SERRATED HANDLE. SMOOTH CURVED SHAFT. (CL-FR-15)



CHALAZION FORCEPS : 95 mm OVERALL LENGTH 8 mm INSIDE DIAMETER OF OPEN UPPER PLATE 10 mm DIAMETER OF SOLID LOWER PLATE SERRATED HANDLE WITH LOCKING THUMB SCREW. (CL-FR-16 /a)



CHALAZION FORCEPS : 95 mm OVERALL LENGTH 15 mm INSIDE DIAMETER OF OPEN UPPER PLATE 17 mm DIAMETER OF SOLID LOWER PLATE SERRATED HANDLE WITH LOCKING THUMB SCREW. (CL-FR-16 /b)



DESMARRES CHALAZION FORCEPS : 95 mm. OVERALL LENGTH 21 X 14 mm. INSIDE DIMENSION OF UPPER PLATE. 24 mm. PLATE WIDE LOWER PLATE. SERRATED HANDLE WITH LOCKING THUMB SCREW. (CL-FR-16 /c)



CASTROVIEJO CORNEAL SCISSORS: CURVED : OVERALL LENGTH 113 mm. SHARP TIP FROM MID SCREW TO TIP 13 mm BLADE. FLAT THREE HOLE BODY. (CL-SC-02)




25

CASTROVIEJO CORNEAL SCISSORS: GENTLY CURVED: OVERALL LENGTH 113 mm. SMOOTH TIP. FROM MID SCREW TO TIP 13 mm BLADE. FLAT THREE HOLE BODY. (CL-SC-02)




26

CASTROVIEJO CORNEAL SCISSORS: STRAIGHT: OVERALL LENGTH 113 mm. FROM MID SCREW TO TIP 13 mm BLADE. SHARP TIP. FLAT THREE HOLE BODY. (CL-SC-02/ST)




27

CASTROVIEJO CORNEAL SCISSORS. ANGLED: OVERALL LENGTH 113 mm. SHARP TIP. 11 mm LONG ANGLED BLADES. FLAT THREE HOLE BODY. (CL-SC-02 / AN)



28

CASTROVIEJO CORNEAL SCISSORS. Right: 7 mm. LONG BLADE. SERRATED / THREE HOLE BODY. (CL-SC-04 /RS) & (CL-SC-04 /RH).




CASTROVIEJO CORNEAL SCISSORS. *Left* : 7 mm. LONG BLADE. SERRATED / THREE HOLE BODY. (CL-SC-04 /LS) & (CL-SC-04 /LH).



WESTCOTT SCISSORS. *Serrated handle* : OVEALL LENGTH 110 mm. BLUNT TIP. FROM MID SCREW TO TIP 14 mm CURVED BLADE. (CL-SC-05 /S-ST)



WESTCOTT SCISSORS. *Perforated handle* : OVERALL LENGTH 110 mm. BLUNT TIP. FROM MID SCREW TO TIP 14 mm CURVED BLADE. (CL-FR-05 /P)



WESTCOTT TENOTOMY SCISSORS : GENTLY CURVED 21 mm BLADE. SERRATED HANDLE. (CL-SC-05-S /CU).



VANNAS SCISSORS, STRAIGHT: OVERALL LENGTH 89 mm. SHARPE TIP.
FROM MID SCREW TO TIP 9 mm BLADE. FLAT SERRATED HANDLE. (CL-SC-03 /SH).




VANNAS SCISSORS: GENTLY CURVED OVERALL LENGTH 89 mm.
FROM MID SCREW TO TIP 11 mm BLADE. SHARP TIP. FLAT SERRATED HANDLE. (CL-SC-03).



ANTERIOR CAPSULOTOMY, TAG CUTTING SCISSORS: OVERALL LENGTH 87 mm.
SPATULATE BLADE. FROM MID SCREW TO TIP 12 mm BLADE. SHARP TIP.
FLAT SERRATED HANDLE. (CL-SC-06).




IRIS SCISSORS: STRAIGHT. OVERALL LENGTH 90 mm. (CL-SC-07).




37

SHARP / BLUNT SCISSORS : STRAIGHT. (CL-SC-08).



38

CASTROVIEJO NEEDLE HOLDER: OVERALL LENGTH -115 mm. FINE 15 mm JAW. POLISHED FINISH. FLAT SERRATED BODY. (CL-NH-01)



39

BARRAQUER NEEDLE HODER. Micro: OVERALL LENGTH -115 mm. FINE 11 mm JAW. POLISHED FINISH. ROUND SERRATED BODY. (CL-NH-02)




40

KALT NEEDLE HOLDER : OVERALL LENGTH 135 mm. FLAT LOCKING HANDLE. HEAVY 12 mm SERRATED STRAIGHT JAWS. (CL-NH-03).



SILCOCK NEEDLE HOLDER : 11 mm JAWS WITH LOCKING HANDLE. (CL-NH-04)



SWISS BLADE BREAKER & HOLDER : 13 mm LONG JAW AUTO SLIDING LOCK.
ROUND KNURLED HANDLE. (CL-BB-01)



KIRBY LENS LOOP : 6 mm WIDE WITH FLAT HANDLE. (CL-MN-01)



SINSKEY LENS MANIPULATOR : 0.2 mm DIAMETER / BLUNT TIP. (CL-MN-02)




CAPSULE POLISHER RING, NIGHTINGALE : 1.75 DIA RING. KNURLING HANDLE. (CL-MN-03)



CAPSULE POLISHER RING, NIGHTINGALE : 1.75 DIA RING. SPIRAL GROOVE HANDLE.
(CL-MN-04)



B.P.HANDLE: BARD-PARKER HANDLE FOR No. 9 KNIFE. KNURLED / FLAT HANDLE.
(CL-BP-01)




FIXATION HOOK, DOUBLE : FLAT / KNURLED HANDLE. 1.5 mm TIP. POLISHED FINISH.
10 mm SPREAD BETWEEN CURVED TIPS. OVERALL LENGTH 5 Inch (Approx.). (CL-MN-05)



KNAPP LACHRYMAL SAC RETRACTOR : BLUNT, 8 mm WIDE FOUR PRONGS. (CL-MN-06)



WELLS ENUCLEATION SPOON : 21 mm. WIDE WITH 5mm BY 8 mm NOTCH. FLAT SMOOTH
HANDLE POLISHED FINISH OVERALL LENGTH 137 mm. (CL-MN-07).



EVISCERATION SPOON : FLAT SMOOTH HANDLE POLISHED FINISH OVERALL LENGTH 137 mm.
8 mm. DIA. (CL-MN-08). 12 mm. DIA. (CL-MN-09).



HEAT CAUTERY COPPER (BALL POINT) : 6 mm COPPER BALL. (CL-MN-10).



HEAT CAUTERY COPPER (TAPER POINT) : TAPER POINT. (CL-FR-10 /T)



BOWMAN LACCHRYMAL PROBE : SIZES "0" to "000". (CL-MN-11)